Airstream' Treated as a construction of the second of the

Class II Biological Safety Cabinets

The World's Most Energy-Efficient, Quiet and Compact Biosafety Cabinet







Removable Paper Catch

Easy to clean

Class II Biological Safety Cabinets

Optional pre-filter can be fitted



Esco Airstream[®] Class II has been certified by PHE / Public Health England (formerly HPA) for compliance to EN 12469

Airflow Sensor

- Monitors real-time airflow for safety
- Alerts the user if airflow is insufficient

Energy-Efficient DC ECM Motor

- The most energy-efficient Class II biosafety cabinet in the world, provides 70% energy savings compared to AC motor
- Stable airflow, despite building voltage fluctuations & filter loading
- Night Setback mode to further reduce power consumption by 60%



З

Power Consumption in Watts for 1.2 m (4 ft) Cabinets

ESZD

CLASS. WORLDWIDE.



	Biosafety Cabinet	Air Quality	Filtration	Electrical Safety
Standards Compliance	EN 12469, Europe SANS 12469, South Africa	ISO 14644.1 Class 3, Worldwide JIS B9920 Class 3, Japan JIS BS 5295, Class 3, UK	EN-1822 (H14), Europe IEST-RP-CC001.3, USA IEST-RP-CC007, USA IEST-RP-CC034.1, USA	IEC 61010-1, Worldwide EN 61010-1, Europe UL 61010-1, USA CAN / CSA-22.2, No.61010-1

*CFDA certification is exclusive to AB2 models sold in China.



				l Safety Ca			
			AC2-2E8 2010718	AC2-3E8 2010658	AC2-4E8 2010621	AC2-5E8 2010656	AC2-6E8 2010657
Glass Side W	all	230 VAC,			AC2-4G8 2010734		AC2-6G8 2010743
		50/60 Hz	AC2-2S8 2010767	AC2-3S8 2010721	AC2-458 2010711	AC2-558 2010725	AC2-6S8 2010722
Cabinet Stainless Ste	el Side Wall				AC2-4D8 2010733		AC2-6D8 2010742
Glass Side W	all	115 VAC, 50/60 Hz	AC2-2E9 2010777	AC2-3E9 2010779	AC2-4E9 2010697	AC2-5E9 2010784	AC2-6E9 2010787
Stainless Ste	el Side Wall		AC2-2S9 2010790	AC2-3S9 2010792	AC2-4S9 2010744	AC2-5S9 2010797	AC2-6S9 2010800
Anti-blowba	Anti-blowback Valve 10"				ABBV-10P 5170352		
Exhaust Ducting Tri-safe Exha	Tri-safe Exhaust Collar with Alarm		N/A		TEM-4 2010606	NIZA	
Thimble Exh	Thimble Exhaust Collar		ECO-AC22 5170520	ECO-AC23 5170521	ECO-AC24 5170522	ECO-AC25 5170523	ECO-AC26 5170524
UV Lamp			UV-1 5170	5A-L)251		UV-30A-L 5170255	
IV Bar	IV Bar		IV-605 5170498	IV-910 5170499	IV-1215 5170231	IV-1520 5170500	IV-1825 5170501
Vork Zone Multiple Piec	e Tray Option	(AC2-S / AC2-D)	SDT-AC2-2E 5020643	SDT-AC2-3E 5020635	SDT-AC2-4E 5020606	SDT-AC2-5E 5020640	SDT-AC2-6E 5020592
Single Piece	Tray Option (A	C2-E / AC2-G)	SGT-AC2-2S 5020696	SGT-AC2-3S 5020648	SGT-AC2-4S 5020627	SGT-AC2-5S 5020651	SGT-AC2-6S 5020645
Pre-filter	Pre-filter		PF-40 5090060	PF-41 5090061	PF-42 5090062	PF-43 5090063	PF-44 5090064
ectrical utlet Direct Moun	Direct Mounted / GFCl				EO		

Service	EU SF-Gas-20 mm		S	F-1G20 517041	0			
	EU SF-Vacuum-20 mm		SF-1V20 5170457					
	EU SF-Air-20 mm		SF-1A20 5170502					
Fixtures	EU SF-Nitrogen-20 mm		SF-1N20 5170503					
	EU SF-Water-20 mm		SF-1W20 5170458					
-	EU SF-Universal-22 mm		SF-2U22 5170504					
Support Stands, Ships Flat	Fixed Stand with Leveling Feet, 28" height	SAL-2A0 Gen 2 5130169	SAL-3A0 Gen 2 5130170	SAL-4A0 Gen 2 5130134	SAL-5A0 Gen 2 5130171	SAL-6A0 Gen 2 5130172		
	Fixed Stand with Leveling Feet, 34" height	SAL-2B0 Gen 2 5130173	SAL-3B0 Gen 2 5130174	SAL-4B0 Gen 2 5130175	SAL-5B0 Gen 2 5130176	SAL-6B0 Gen 2 5130177	ABBV-10P	
	Fixed Stand with Caster Wheels, 28" height	SPC-2A0 Gen 2 5130161	SPC-3A0 Gen 2 5130155	SPC-4A0 Gen 2 5130152	SPC-5A0 Gen 2 5130162	SPC-6A0 Gen 2 5130154		
	Fixed Stand with Caster Wheels, 34" height	SPC-2B0 Gen 2 5130164	SPC-3B0 Gen 2 5130165	SPC-4B0 Gen 2 5130166	SPC-5B0 Gen 2 5130167	SPC-6B0 Gen 2 5130168		
	Telescopic Stand with Leveling Feet, 1" adjustment	STL-2A0 5130092	STL-3A0 5130050	STL-4A0 5130051	STL-5A0 5130052	STL-6A0 5130053	SF	
	Telescopic Stand with Caster Wheels, 1" adjustment	STC-2A0 5130135	STC-3A0 5130055	STC-4A0 5130056	STC-5A0 5130057	STC-6A0 5130058		
	Motorized Height Stand with Caster Wheels, 39.5" height		SPM-3A2 5130093	SPM-4A2 5130047	SPM-5A2 5130100	SPM-6A2 5131141	-	
Misc	IQ/OQ Protocol		9010179					

AC2 Biological Safety Cabinet Engineering Drawing



- 4. Stainless Steel Work Tray (available in single and multi pieces)
- 5. Exhaust Collar (optional)
- 6. Exhaust ULPA / H14 Filter

- 10. T5 Fluorescent Lamps
- 11. IV Bar Retrofit Kit Provision
- 12. Electrical outlet Retrofit Kit Provision (1 on each side)
- 13. Service Fixture Retrofit Kit Provision (2 on each side)

Optional Exhaust Collar Positions for Thimble-Ducting for AC2 Models



C :	2	3	4	5	6	ft	
Size	0.6	0.9	1.2	1.5	1.8	m	
х	233	331	408	560	560		
Y	334	334	334	334	326	mm	
х	9.2	13	16	22	22	inches	
Y	13.1	13.1	13.1	13.1	12.8	inches	

18. UV Lamp (optional)

19. Removable Paper Catch (with optional pre-filter)

SPC-_A0 Gen2

5

Glass Side: 230 VAC, 50	0/60 Hz	AC2-2E8 2010718	AC2-3E8 2010658	AC2-4E8 2010621	AC2-5E8 2010656	AC2-6E8 2010657		
Glass Side: 115 VAC, 50/60 Hz		AC2-2E9 2010777	AC2-3E9 2010779	AC2-4E9 2010697	AC2-5E9 2010784	AC2-6E9 2010787		
Stainless Steel Side: 23	80 VAC, 50/60 Hz	AC2-258 2010767	AC2-358 2010721	AC2-458 2010711	AC2-558 2010725	AC2-6S8 2010722		
Stainless Steel Side: 115 VAC, 50/60 Hz		AC2-259 2010790	AC2-359 2010792	AC2-459 2010744	AC2-559 2010797	AC2-659 2010800		
Nominal Size		2 ft (0.6 meter)	3 ft (0.9 meter)	4 ft (1.2 meter)	5 ft (1.5 meter)	6 ft (1.8 meter)		
	Width	730 mm (28.8")	1035 mm (40.8")	1340 mm (52.8")	1645 mm (64.8")	1950 mm (76.8")		
External Dimensions	Depth without Arm Rest		-	753 mm (29.5")	-			
(W x D x H)	Depth with Arm Rest	810 mm (32.0")						
	Height	1400 mm (54.8")						
C	Width	610 mm (24.0")	915 mm (36.0")	1220 mm (48.0")	1525 mm (60.0")	1830 mm (72.0")		
Gross Internal Dimensions	Depth	580 mm (22.8")						
(W x D x H)	Height	660 mm (26.0 °)						
Usable Work Area		0.27 m ² (2.9 sq.ft.)	0.42 m ² (4.5 sq.ft.)	0.56 m² (6.1 sq.ft.)	0.71 m ² (7.63 sq.ft.)	0.86 m ² (9.2 sq.ft.)		
Tested Opening			1	175 mm (7 ")	-			
Work Opening		 190 mm (7.5")						
Average Airflow	Inflow	0.45 m/s (90 fpm)						
Velocity	Downflow	0.30 m/s (60 fpm)						
	Inflow	173 cmh (102 cfm)	259 cmh (152 cfm)	346 cmh (204 cfm)	432 cmh (254 cfm)	519cmh (305 cfm)		
	Downflow	369 cmh (217 cfm)	553 cmh (325 cfm)	738 cmh (434 cfm)	922 cmh (543 cfm)	1107 cmh (657 cfm)		
Airflow Volume	Exhaust	173 cmh (102 cfm)	259 cmh (152 cfm)	346 cmh (204 cfm)	432 cmh (254 cfm)	519cmh (305 cfm)		
Airnow volume	Required Exhaust with Optional Thimble Exhaust Collar	260 m³/h (153 cfm)	320 m³/h (189 cfm)	538 m³/h (317 cfm)	615 m³/h (362 cfm)	823 m³/h (485 cfm)		
	Static Pressure for Optional Thimble Exhaust Collar	28 Pa / 0.11 in H ₂ O	29 Pa / 0.11 in H ₂ O	31 Pa / 0.12 in H ₂ O	35 Pa / 0.14 in H ₂ O	47 Pa / 0.18 in H ₂ O		
ULPA Filter Typical Effic	iency	>99.999% at 0.1 to 0.3 micron, ULPA as per IEST-RP-CC001.3 USA						
		T	>99.999	% at MPPS, H14 as per EN	N 1822 EU			
Sound Emission*	NSF / ANSI 49	56.3	56.6	58.7	58.2	59.4		
	EN 12469	51.0	52.0	53.5	53.6	55.7		
Fluorescent Lamp Inten	sity (Lux)	859	1279	1404	1227	1384		
Fluorescent Lamp Inten	sity (foot-candles)	80	119	130	114	129		
	Main body	1.2 mm (0.05") 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester Isocide™ antimicrobial powder-coated finish						
Cabinet Construction	Work Zone		1.5 mm (0.06") 16 g	gauge stainless steel, typ	e 304, with 4B finish			
	Side Walls (E Series)		UV-absorbing tempere	ed glass, 5 mm (0.2"), co	lorless and transparent			
	Side Walls (S Series)		1.5 mm (0.06") 16 g	gauge stainless steel, typ	e 304, with 4B finish			
Electrical	Cabinet Full Load Amps (FLA)	1.8	3.5	3.7	4.3	5.5		
	Heat Load (BTU/Hr)	324	447	580	717	966		
Nominal Power Consun	nption (W)	95	131	160	210	283		
Net Weight**		116 Kg (256 lbs)	173 Kg (381 lbs)	230 Kg (507 lbs)	288 Kg (635 lbs)	346 Kg (763 lbs)		
Shipping Weight**		143 Kg (315 lbs)	214 Kg (472 lbs)	285 Kg (628 lbs)	356 Kg (785 lbs)	428 Kg (944 lbs)		
Shipping Dimensions Maximum (W x D x H)*	*	850 x 820 x 1760 mm (33.5" x 32.3" x 69.3")	1120 x 820 x 1760 mm (44.1" x 32.3" x 69.3")	1450 x 820 x 1760 mm (57.1" x 32.3" x 69.3")	1720 x 820 x 1760 mm (67.7" x 32.3" x 69.3")	2050 x 820 x 1760 mn (80.7" x 32.3" x 69.3"		
Shipping Volume, Maximum**		1.23 m ³ (43.2 ft ³)	1.62 m³ (57.2 ft³)	2.09 m ³ (73.8 ft ³)	2.48 m³ (87.6 ft³)	2.96 m ³ (104.5 ft ³)		

*Noise reading in open field condition / **anechoic** chamber. Noise reading in **normal room varies** by room size, layout, and background noise, but may reach roughly 3-4 dBA above these values. **Cabinet only, excludes optional stand.

TECHNICAL SPECIFICATIONS						
Glass Side: 230 VAC, 50	1/60 Hz	AC2-4G8 2010734	AC2-6G8 2010743			
Stainless Steel Side: 23	0 VAC, 50/60 Hz	AC2-4D8 2010733	AC2-6D8 2010742			
Nominal Size		4 ft (1.2 meter)	6 ft (1.8 meter)			
	Width	1340 mm (52.8")	1950 mm (76.8")			
External Dimensions	Depth without Arm Rest	753 mm (29.5")				
(W x D x H)	Depth with Arm Rest	810 mm	(32.0")			
	Height	1400 mm (54.8")				
Gross Internal	Width	1220 mm (48")	1830 mm (72")			
Dimensions	Depth	580 mm (22.8")				
(W x D x H)	Height	660 mm (26")				
Usable Work Area		0.56 m² (6.1 sq.ft.)	0.86 m² (9.0 sq.ft.)			
Tested Opening		175 mi	m (7 ")			
Working Opening		190 mm	ו (7.5")			
Average Airflow	Inflow	0.45 m/s	(90 fpm)			
Velocity	Downflow	0.30 m/s (60 fpm)				
	Inflow	346 cmh (588 cfm)	519 cmh (881 cfm)			
	Downflow	738 cmh (1254 cfm)	1107 cmh (1880 cfm)			
	Exhaust	346 cmh (588 cfm)	519 cmh (881 cfm)			
Airflow Volume	Required Exhaust with Optional Thimble Exhaust Collar	538 m³/ h (317 cfm)	823 m³/ h (485 cfm)			
	Static Pressure for Optional Thimble Exhaust Collar	31 Pa / 0.12 in H ₂ O	47 Pa / 0.18 in H ₂ O			
		>99.999% at 0.1 to 0.3 micron, L	JLPA as per IEST-RP-CC001.3 USA			
ULPA Filter Typical Effici	ency	>99.999% at MPPS, H14 as per EN 1822 EU				
	NSF / ANSI 49	61.3 dBA	62.5 dBA			
Sound Emission*	EN 12469	58.3 dBA	59.5 dBA			
Fluorescent Lamp Intens	sity (Lux)	1400				
Fluorescent Lamp Intens	sity (foot-candles)	130				
	Main body	1.2 mm (0.05") 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester Isocide™ antimicrobial powder-coated finish				
Cabinet Construction	Work Zone	1.5 mm (0.06 ") 16 gauge stainless steel, type 304, with 4B finish				
	Side Walls (G-Series)	UV-absorbing tempered glass, 5 mi	m (0.2 "), colorless and transparent			
	Side Walls (D-Series)	1.5 mm (0.06") 16 gauge stainle	ss steel, type 304, with 4B finish			
	Cabinet Full Load Amps (FLA)	9.4 A	12.6 A			
Electrical	Heat Load (BTU / Hr)	905	1230			
Nominal Power Consumption		265 W	360 W			
Net Weight**		240 Kg (529 lbs)	366 Kg (807 lbs)			
Shipping Weight**		295 Kg (650 lbs)	448 Kg (988 lbs)			
Shipping Dimensions, N	laximum (W x D x H)**	1450 x 820 x 1760 mm (57.1" x 32.3" x 69.3")	2050 x 820 x 1760 mm (80.7" x 32.3" x 69.3")			
Shipping Volume, Maxir	mum**	2.09 m³(73.8 ft³)	2.96 m ³ (104.5 ft ³)			

*Noise reading in open field condition / **anechoic** chamber. Noise reading in **normal room varies** by room size, layout, and background noise, but may reach roughly 3-4 dBA above these values. **Cabinet only, excludes optional stand.

First Airstream[®] Offers the Most Complete Class II Cabinet Range

Airstream Product	E-Series	G-Series	S-Series	D-Series		
Side Wall Tempered glass increases visibility and prevents the operator from experiencing a "boxed-in" feeling			One-piece stainless steel with coved corners for cleanability. Side capture zones and negative pressure side walls optimize containment.			
Work Tray	Multi-piece,	Autoclavable	Single-piece stainles	s steel, spill retaining		
Fan System	Single blower for inflow and downflow. Energy-efficient and cost-effective	Dual blowers for inflow and downflow. Redundant system provides protection in case of fan failure.	Single blower for inflow and downflow. Energy-efficient and cost-effective.	Dual blowers for inflow and downflow. Redundant system provides protection in case of fan failure.		
Exhaust Filter	Cost-effective ULPA/H14 filter with > 99.999% efficiency	Dual ULPA/H14 filters that provide > 100x better protection than single filter system	Single ULPA/H14 Filter >99.999% Efficient, Cost Effective	Dual ULPA/H14 Filters, >100.000x Better Protection than Single Filter System		
Size Available	0.6 m (2'), 0.9 m (3'), 1.2 m (4'), 1.5 m (5'), 1.8 m (6')	1.2 m (4'), 1.8m (6')	0.6 m (2'),0.9 m (3'), 1.2 m (4'), 1.5 m (5'), 1.8 m (6')	1.2 m (4'), 1.8 m (6')		

AC2-D Airflow Diagram



ULPA-filtered air

Unfiltered / potentially contaminated air

Room air / Inflow air

AC2-G Airflow Diagram



ULPA-filtered air

Unfiltered / potentially contaminated air

Room air / Inflow air

Comprehensive Performance Testing At Esco

Every Airstream® AC2 model manufactured by Esco is individually tested, documented by serial number and validated with the following test methods:

- Inflow and downflow velocity
- PAO aerosol challenge for filter integrity
- Airflow pattern visualization
- Electrical safety to IEC61010-1

 Additional KI-Discus containment and microbiological testing are performed on statistical sampling basis.





ESCO GLOBAL NETWORK 43 LOCATIONS IN 23 COUNTRIES ALL OVER THE WORLD



Regional Distribution Centers



ART Equipment **Biological Safety Cabinets** CO₂ Incubators Compounding Pharmacy Equipment Containment / Pharma Products Ductless Fume Hoods Lab Animal Research Products Laboratory Centrifuges Laboratory Fume Hoods Laboratory Ovens and Incubators Laboratory Shakers Laminar Flow Clean Benches PCR Cabinets PCR Thermal Cyclers Powder Weighing Balance Enclosures Ultra-low Temperature Freezers



escoglobal.com

ISOCIDE[™]

The Esco Group of Companies is a global life sciences tools provider with sales in over 100 countries. The group is active in lab equipment, pharma equipment and medical devices. Manufacturing facilities are located in Asia and Europe. R&D is conducted worldwide spanning the US, Europe and Asia. Sales, service and marketing subsidiaries are located in 12 major markets including the US, UK, Singapore, Japan, China and India. Regional distribution centers are located in the US, UK, and Singapore.

Life Science • Chemical Research • Assisted Reproductive Technology (ART) • Pharmaceutical Equipment • General Equipment



WORLD CLASS. WORLDWIDE.

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com www.escoglobal.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA Tel 215-441-9661 • Fax 484-698-7757 eti.admin@escoglobal.com • www.escolifesciences.us

Esco Global Offices: Bangladesh | Cameroon | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | United Kingdom | USA | Vietnam





